
SARS-CoV-2 Nucleocapsid (N) Antigen Rapid Detection Kit


(Colloidal gold method)

User's Manual


For Self-test Use Only! Please read this manual carefully before operation to ensure proper use.

PRODUCT NAME:SARS-CoV-2 Nucleocapsid (N) Antigen Rapid Detection Kit
(Colloidal gold method)**FOR AGES:** 15-70**PACKING SPECIFICATIONS:** 1/2/3/5/6/7/8/9/10/12/15/20/25/30/50 test (s) /kit**INTENDED USE:**The product is used for the qualitative detection of SARS-CoV-2 infection. SARS-CoV-2 virus can cause death and complications. No instrument required. It is used for the screening of infected patients. This kit is used for personal use by untrained laymen.**TEST PRINCIPLE:**This product uses colloidal gold immunochromatography combined with the double-sandwich principle to detect the N-antigen of novel coronavirus in human nasal swab samples. After adding the sample to be tested, the SARS-CoV-2 antigen in the sample and the coated colloidal gold-IgY-SARS-CoV-2 antibody-1 conjugate on the binding pad undergo specific immunological binding to form an antigen-colloidal gold complex. As the liquid flows, the antigen-colloidal gold complex binds to the novel coronavirus antibody-2 coated on the detection area (T) of the nitrocellulose membrane, and is captured by the detection area (T), forming a color band. The IgY colloidal gold marker is specifically immunologically bound to the anti-IgY antibody coated on the nitrocellulose membrane quality control area (C). The quality control area (C) is colored whether or not the sample contains novel coronavirus antigens.**STORAGE CONDITIONS AND VALIDITY PERIOD:**


- Storage conditions: The original packaging should be stored in a dry place at 4~30°C, protected from light, and do not freeze.
- Validity period: 18 months.
- The test pad should be used as soon as possible within 1 hour after unpacking the aluminum foil bag.

TEST METHOD:**INTERPRETATION OF TEST RESULTS:****Positive:**

Two red bands, the detection area (T line) and the quality control area (C line) are both colored.

Note:A positive result means it is likely that you have COVID-19. To confirm the result, repeating the test after 24 to 36 hours is recommended. Test results should always be considered in the context of clinical observations and epidemiological data in making final diagnostic and patient management decisions. If you are tested positive, immediately self-quarantine following your local guidelines and contact your general practitioner/doctor or the local health department in accordance with the instructions of your local authorities.**METHODE DE TEST:****INTERPRETATION DES RESULTATS DE TEST****Positif**

Two red bands, the detection area (T line) and the quality control area (C line) are both colored.

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2. Interference:

The following internal and external interfering substances were tested, and the test results were negative, indicating that there is no interference with the following test substances.

| Analytes | Concentration | Analytes | Concentration |
|-----------------|---------------|-----------------------------|---------------|
| Purified murine | s1.0mg/L | Sodium chloride | s0.9% |
| α-interferon | s100U/mL | Fluorhexidine | s200µg/mL |
| Rituximab | s2.0mg/L | Budesonide | s100µM/L |
| Oscilatumab | s375µg/L | Benzocaine | s150mg/L |
| Perimavimab | s20µg/mL | Mupirocin | s1.05mg/L |
| Levoferrocim | s0.5µg/mL | Dihydroxyacetone | s0.075mg/L |
| Cefazosrine | s100mg/mL | Devotoximab | s0.0105mg/L |
| Meropenem | s1.0µg/mL | Desmethylmethane | s1.2mg/L |
| Tobramycin | s0.125mg/mL | Normal nasopharyngeal mucus | s5% |
| Oxymecetazole | s0.125µg/mL | Ganze Blut | s5% |

3. Clinical performance

Compare the results of SARS-CoV-2 Nucleocapsid (N) Antigen Rapid Detection Kit with the results of novel coronavirus nucleic acid detection reagent (RT-PCR) in valid cases. The results are shown below.

Statistics of test results of nasal swab samples:

| Reagents | RT-PCR | Positive | Negative | Total |
|--------------------------------------------------------------------------------------|--------|----------|----------|-------|
| SARS-CoV-2 Nucleocapsid(N) Antigen Rapid Detection Kit (Colloidal Gold Method) | | 105 | 0 | 105 |
| | | 5 | 456 | 461 |

| Statistical Indicators | Estimated value | 95%CI |
|---------------------------|-----------------|-------------------|
| Positive coincidence rate | 95.45% | 89.71% ~ 98.51%</ |

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