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Real Time PCR • LAMP PCR • RAPID Ab •
RAPID Ag • RAPID Ag FIA

GENEDIA COVID-19 Ag

GENEDIA W COVID-19 Ag device is a chromatographic immunoassay for the qualitative detection of specific antigens to COVID-19 present in human nasopharyngeal swab and sputum.

Immediate on-site Antigen testing

- Allow wider testing with fast test time (10 minutes)
- All necessary materials provided / no equipment needed
- High capacity to meet the most urgent medical and public health needs

Performance compared with other methods

For COVID-19 Ag		Comparator (RT-PCR)		Total
		Positive (+)	Negative (-)	
GENEDIA W COVID-19 Ag	Positive (+)	182	0	182
	Negative (-)	20	229	249
Total		202	229	431




Positive results broken down by days since symptom onset

Days Since Symptom Onset	RT-PCR Positive (+)	GENEDIA W COVID-19 Ag Positive (+)	PPA (Positive Predictive Agreement)
≤7	68	64	94%
Asymptomatic	54	46	85%

- Clinical Sensitivity = 90.1% (95% CI : 85.2% - 93.5%)
- Clinical Specificity = 100.0% (95% CI : 98.4% - 100%)
- Agreement : 95.3%

Made in Korea

TEST PROCEDURE

				
1 STEP Open the extraction solution tube.	2 STEP Insert the patient swab sample into the extraction solution tube. Then, swirl swab at least 6 times.	3 STEP Roll the swab head against the inside of the tube as you remove it.	4 STEP Place the sample developing screw cap into the extraction tube and close it.	5 STEP Add 3 drops of the extraction solution with specimen into the sample well. After 10 minutes, interpret the test results.

INTERPRETATION OF THE RESULTS

COVID-19 Ag	COVID-19 Ag	COVID-19 Ag	COVID-19 Ag
			
Negative	Positive	Invalid	


Positive Result

Ag	Ag	Ag	Ag	Ag
				

If there is any blurry line or color change on the T line within 15 minutes, it should be read as positive.

SPECIMEN COLLECTION

Nasopharyngeal Swab

Sputum


1 Rinse your mouth with water

2 Breath in and out 3 times

3 Give a sputum sample

4

Put the swab into the container containing the sputum, turn it several times and soak it sufficiently

MATERIALS PROVIDED


SPECIFICATION

Method	Immunochromatography
Packing unit	20 Tests / kit
Certification	CE-IVD
Specimen	Nasopharyngeal swab and sputum
LOD	7.50×10^2 TCID ₅₀ /mL
Running time	10 minutes
Expiry date	12 months from the date of manufacture
Storage condition	2 ~ 30°C

GENEDIA  Series
RAPID Ag

- ✓ Specimen by Nasopharyngeal Swab or **SPUTUM**
- ✓ Reading Time only **5~10min**
- ✓ LOD(Limit of Detection) **7.50×10^2 TCID₅₀/mL**
- ✓ Reliable Manufacturer  **Made in KOREA**
- ✓ Germany **BfArM Listed**  Bundesinstitut für Arzneimittel und Medizinprodukte

Sample material for direct pathogen detection *by RKI at 25.Jan.2021*

RKI (Robert Koch Institute, Germany) comments the proper example of Sample material for direct pathogen detection as "Nasopharyngeal swab", "Oropharyngeal swab (throat swab)", "Bronchoalveolar lavage", "**Sputum**(in patients with a productive cough)".
The use of other sample materials, such as throat rinse water / gargle water and saliva, is discussed on various occasions.

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WHO IS GCMS(Green Cross Medical Science Corp.)?