DrugSIP[®]

Multi Analyzer

Drug detection device

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CONTENTS

1- DEVICE APPLICATIONS.....	6
2- DESCRIPTION.....	6
3- TECHNICAL SPECIFICATIONS.....	10
4- UNPACKING AND PREPARING.....	11
5- INSTRUCTIONS FOR USE.....	14
6- TEMPERATURE CONTROL OF SAMPLE.....	30
7- SOFTWARE PC.....	31
8- TROUBLE SHOOTING.....	36
9- READING DEVELOPED TESTS.....	37
10- SAFETY.....	39
11- MAINTENANCE AND STORAGE.....	40C
12- LEANING.....	41

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DrugSIP®

Multi Analyzer

The team at Alpha San Ignacio Pharma (AlphaSIP) Laboratories would like to thank you for purchasing our DrugSIP® Multi Analyzer, designed for the detection of multiple drugs in saliva and urine. We hope you like the product and find it useful. We recommend you care for the equipment according to the guidelines in this manual.

AlphaSIP is a Spanish company pioneering in the development of medical diagnosis sensors which can be used in various fields, such as road safety and occupational health. We have used the most advanced electronic technology for rapid detection of recreational drugs in saliva and urine, with the most hygienic and reliable sample collection, using our DrugSIP® Multi Analyzer device. AlphaSIP products are developed according to the EC marking guidelines, and are based on a robust ergonomic design and a commitment to protecting the user's safety. The quality of the materials used and the proper use of the device will prolong its useful life. Incorrect or inappropriate use of the equipment may result in errors or damage. Please read the Maintenance points (Point 11), which include safety aspects. AlphaSIP reserves the right to make improvements and/or changes to the products and documentation there of, in order to update, revise or modify the products and/or accessories.

PLEASE READ THIS INSTRUCTION MANUAL CAREFULLY BEFORE USING THE EQUIPMENT, IN ORDER TO MAXIMISE ITS PERFORMANCE AND PROLONG ITS LIFE

Be particularly aware of the following:

- This manual is an inseparable part of the DrugSIP® Multi Analyzer equipment, and therefore must be available to all users of the equipment.
- The device must always be handled carefully, avoiding striking it, dropping heavy or sharp objects on it or dropping the device itself, as well as taking care not to spill liquids on it.
- Never dismantle the equipment in order to mend it yourself; in addition to invalidating the guarantee, this could also produce malfunction of all the equipment, as well as harming the people handling it.

- Your distributor can clarify any doubts (start-up, operation). You can also email your questions or suggestions to us at AlphaSIP Technical Service: contact@alphasip.es
- This equipment is covered by the Royal Legislative Decree 1/2007, of 16 November, of the General Law for the Protection of Consumers and Users.
- Equipment inspections which take place after the first year of purchase are not included in the guarantee.
- Handling of the equipment by unauthorised staff will result in total loss of the guarantee cover.
- Accessories and replacements, as well as their loss, are not covered by the guarantee. Normal wear and tear of parts is also not covered in the guarantee period.
- Ensure you keep the purchase invoice in case you need to make a complaint or make use of the guarantee. If you send the equipment to Technical Service enclose the invoice or a copy of it as a guarantee document.
- The manufacturers reserve the right to make changes and improvements to this manual and equipment.

WARNING!

ALPHASIP ACCEPTS NO LIABILITY FOR DATA OBTAINED FROM INCORRECT HANDLING OF THE DEVICE

1. Device applications

The DrugSIP® Multi Analyzer allows fast and effective detection of the presence of recreational drugs in saliva and urine. The system carries out the whole process automatically and identifies the presence of drugs by reading the immune-reactive strips in the DrugSIP® Cassette. It is intended for use by the security forces, and in laboratories, pharmacies, schools, hospitals, doctors' surgeries and road safety.

2. Description

The test to detect drug use consists of two main parts. Collecting the sample using the DrugSIP® Cassette (saliva or urine), and a second stage, automatic processing and analysis of the results using the DrugSIP® Multi Analyzer electronic processing and reading device.

2.1 Oral Fluid DrugSIP® Cassette kit

The Oral Fluid DrugSIP® Cassette kit provided by AlphaSIP consists of the following parts:

- Swab (cotton bud). Appliance with a saturation indicator for collecting and storing the optimal amount of saliva sample for performing the test. It should be found in an individual bag.
- Drug detection kit. This includes the diagnosis card/cassette containing the immune-reactive test strips and the base that hold the conjugated antibodies with colloidal gold and a plastic vial with the liquid (buffer) for the detection. It should be found in an individual bag.

Figure 1 illustrates the components of the saliva sample collection kit Oral Fluid DrugSIP® Cassette:


Figure 1. Components of the Oral Fluid DrugSIP® Cassette sample collection kit

WARNING!

Do not touch the card with the test strips. It should only be touched if its original position has changed before the test. In this case, it must be repositioned correctly, holding the card at the plastic edges. Before using, read the instructions included with the Oral Fluid DrugSIP® Cassette.

2.2 Urine DrugSIP® Cassette kit

The Urine DrugSIP® Cassette kit provided by AlphaSIP consists of the following items:

DrugSIP® urine collection kit. It is included in an individual bag and it consists of:

- Sterile urine collection container.
- Thermometric strip to verify the reliability of the obtained sample.
- Protection sticker to proof the sample contact area

DrugSIP® urine detection kit. It is included in an individual bag and it consists of:

- Diagnosis card/cassette containing the immuno-reactivate strips.
- Base for placing the diagnosis card to read the test automatically.
- Plastic pipette for dispensing the sample.

Figure 2 shows the components of the Urine DrugSIP® Cassette kit:

Figure 2. Components of the Urine DrugSIP® Cassette sample collection kit

WARNING!

Take care not to touch the reactive strips under any circumstances. When you are handling the card for dispensing the urine sample or placing it in the base, avoid contact. Carry out these procedures by holding the card by its plastic edges. Carefully read the instructions included with the Urine DrugSIP® Cassette before use.

2.3 DrugSIP® Multi Analyzer device

The DrugSIP® Multi Analyzer (Figure 4) is a kit for detecting drugs of abuse. It comes in a carrying case for protection and easy, convenient, safe transport.

The DrugSIP® Multi Analyzer allows to process and read the test (Blue button) or only to read the test (Yellow button).

The contents of the carrying case can be seen in Figure 3 and it consists of the following components:

- DrugSIP® Multi Analyzer (electronic processing, reading and documentation device).
- Carrying case.
- External power source.
- Mini-USB Cable for data access.
- Instruction manual.
- 2 rolls of paper for thermal printer.
- Software installation CD for data downloading.


Figure 3. Contents of the carrying case


Figure 4 DrugSIP® Multi Analyzer

3. Technical specifications

DrugSIP® Multi Analyzer electronic device	
Optical reader	Yes
Reading time	5-7 min
Full charge time	6-8 hrs
Full charge runtime	6 hrs
Battery status indicator	LED display and on home screen
Test memory	Last 10.000 tests
Temperature control of sample	Yes
LCD Display	2 rows with 16 characters
Battery charger mains voltage	100-240 V AC
Battery charger mains frequency	50/60 Hz
Battery charger output voltage	12 V DC
Battery charger output current	2A
Dimensions (Height x Width x Depth)	358.5 x 180 x 195 mm
Weight	2,3 Kg
Integrated printer	Yes
Test printing time	14 s
Full charge runtime	150 prints-out

4. Unpacking and Preparing

4.1 Unpacking and checking material

After received the goods provided by AlphaSIP , verify that the package has not been manipulated or beaten and proceed as follows:

- Unpack the equipment.
- Without switching on the equipment, ensure that no damage has occurred in transit. If so, contact your distributor immediately in order to make a complaint within the time period established by the transport service.
- Check that you have received all the components included in the case.
- The equipment can be returned within the 15 days following dispatch, provided it arrives complete, in its original packaging and with all accessories and documents included.

4.2 Preparing the device

Before using the device it is recommended that you familiarise yourself with its components and basics, as well as how the controls work.

**READ THIS INSTRUCTION MANUAL CAREFULLY BEFORE USING
THE EQUIPMENT, IN ORDER TO MAXIMISE ITS PERFORMANCE AND
PROLONG ITS LIFE**

- Before use, ensure that the battery is sufficiently charged. If in doubt, charge the device for at least 4 hours. For prolonged use of the device we recommend full charging (6-8 hours).

- Once the charge is complete it is advisable to disconnect the device, as it could become damaged.

- Check that the printer has enough paper for the anticipated number of print-outs. Provide yourself with enough replacement paper, supplied by AlphaSIP.

4.2.1 Language setting / Time and Date settings

Before using the DrugSIP® Multi Analyzer the language, the date and the time of the device must be set by performing the following steps:

1. Turn on the DrugSIP® Multi Analyzer while holding down the blue button for at least 5 seconds until the screen shows **"CHANGE LANGUAGE DATE AND TIME"**. (This message will appear in the device's default language)

WARNING: If you have accidentally accessed to this mode, press the yellow button to return to the main screen "DrugSIP Ready for test"


2. To start changing these settings, press the blue button. First, you can choose the language. **"LANGUAGE"** will appear on the main screen. Press the yellow button until the screen shows the language you would like to choose (English / Spanish / Deutsch). Then press the blue button.

3. Later access the date settings. The message **"WEEK DAY"** appears on the screen, then press the yellow button until you see the correct day of the week (Monday, Tuesday...). Then press the blue button to show **"WEEK MONTH"** and press the yellow button until you see the day of the current month. By pressing the blue button again, **"MONTH"** appears on the screen, and as above, use the yellow button to select the month. Press the blue button again and you will see **"YEAR"**. Choose the correct year with the yellow button and press the blue button.

4. Then access the time settings. The message **"HOUR"** will appear on the screen. Use the yellow button to select the current hour. By pressing the blue button again, the screen will show **"MINUTES"**, press the yellow button to change the minutes. By pressing the blue button again, **"SECONDS"** appears on the screen and the seconds setting can be changed with the yellow button.

5. Finally, press the blue button again to perform system boot, the screen will show **"Initiating..."**. A few seconds later the system returns to the main screen **"DrugSIP Ready for test"**.

WARNING: If any button has been pressed and the hour or the date have not been adjusted correctly, the device must be switched off and the previous steps must be performed again.

4.2.2 Replacing printer paper roll

To replace or check the printer paper roller of the DrugSIP® Multi Analyzer, perform the following steps:

1. First, unlock the paper compartment cover by raising the tab from closed position and lift the cover to open it (Figure 5a).
2. Remove the spool from the used paper roll in the paper roll compartment.
3. Loosen the end of the paper, be taped or sealed with an adhesive material. Unroll the loose end 4 or 5 centimetres from the roll.
4. Place the paper roll in the compartment, making sure that the paper feeds correctly from the underside of the roll. This aligns the paper for the thermal print head (Figure 5b).


Figure 5. (a) unlock the paper compartment (b) place the paper roll

5. Close the paper cover, making sure the loose end of the paper roll is outside of the compartment. The tab will come down automatically when the compartment is closed.

WARNING: If the device does not print the ticket, probably the printer paper roll is run out.
WARNING: If after replacing the printer paper roll, the device prints white tickets or without symbols, make sure the paper roll is positioned correctly; as it may be positioned the wrong way round.

5. Instructions for use

The DrugSIP® Multi Analyzer device allows individual tests to be carried out very easily, simply by inserting DrugSIP® Cassette kit in the device and waiting for it to incubate and analyse the test.

WARNING: For correct operation, you must ensure that the device is placed on a stable, flat, horizontal surface without inclination.

5.1 Switching on and battery charging

5.1.1 Switching on

1. First of all, switch on the device using the ON/OFF (Figure 4) switch on the back of the device.

WARNING: When using the device with an external power source, it must be installed near an easily accessible wall socket.

2. The following message will then appear on the screen “DrugSIP ALPHASIP” together with the software version used in the device, eg. “V 1.0”.

3. Then the battery charge level of the device will automatically appear; the message on the screen will be similar to “Battery: 87%”.


4. A few seconds later, the device shows the message “Initiating...” accompanied by a progressive charge bar.

WARNING: The device must not be turned off while the screen shows “Initiating...”, as it could damage the system.

WARNING: Do not perform a test when the battery is going to run out

When the system is operative and the main screen shows “DrugSIP Ready for test”, users can start processing and reading a test.

5. These steps take a few seconds and take place automatically as soon as the device is switched on.

5.1.2 Battery charging

The DrugSIP® Multi Analyzer device is an autonomous system that can be used without being connected permanently to the mains.

Battery indicator leds

RED indicator light

- Turned on. The battery level is low and the device must be connected to a power source.

GREEN indicator light

- Turned on. The device is connected to a power source and the battery is in charging process.
- Turned on and blinking. The device’s batteries are about to be charged at their maximum capacity.
- Turned off and connected to a power source. The device is completely charged.

1. When the “Low battery” message appears on the screen and the red indicator light comes on, this means the battery is going to run out.

In this case, the DrugSIP® Multi Analyzer device must be connected to a power source.

WARNING: Do not perform a test when the red indicator light comes on and the device is not connected to a power source,

WARNING: Do not check the battery level when the device is connected to a power source because this level may not be correct.

2. The device's green light comes on during the battery charging process. Once the charge is complete the green light goes off.

WARNING: Disconnect the DrugSIP® Multi Analyzer device from the mains when it is completely charged. When the device is charging the sound of the internal fan will be heard. This is automatically activated to prevent overheating's device.

INFREQUENT USE OF THE DEVICE

The battery of the DrugSIP® Multi Analyzer needs a frequent maintenance to prevent it from damage.

- Do not store the device for a long period of time without performing a battery maintenance.
- Under no circumstances the battery must be fully discharged.
- If possible the device must be stored with a 100% battery level.

WARNING: The device must be charged at least 3 weeks, provided that the device is always stored with a battery level between 90% and 100%.

5.2 Performing a simple saliva test


5.2.1 Collecting a saliva sample

To collect the saliva sample use the swab included in the Oral Fluid DrugSIP® Cassette provided by AlphaSIP. The kit consists of the saliva collector swab and the detection kit, as seen in Figure 1. To collect the saliva sample follow these steps:


Read the enclosed instructions on Oral Fluid DrugSIP® Cassette before use. They may include improvements not described in this manual.

1. Give the collection swab to the person who is going to take the test. This person should take the swab out of the bag and put it into their mouth. To begin the collection, rub the teeth, gums and all around the mouth (cheeks) with the sponge. They must keep it in their mouth until the sponge is saturated with the saliva sample. To know it is saturated, check that the sponge of the swab has expanded as far as a small slot, taking up almost all the available space. In figure 6 you can see the proper saturation of sample in order to carry out the test properly.

SUFFICIENT SATURATION

INSUFFICIENT SATURATION

Figure 6. Sufficiently saturated sample

2. Open the bag of Oral Fluid DrugSIP® Cassette and take the diagnosis card and the base out of the bag. Check that the buffer vial is filled with liquid (Figure 1). **VERY IMPORTANT:** if the vial does not have liquid inside, dispose the test and use a new one.

3. While the subject collects the sample, his or her identification must be written on the front of the diagnosis card in the space allocated (Figure 7).


Figure 7. Allocated space


Figure 8. Oral Fluid DrugSIP® Cassette (a) swab locating (b) plastic buffer vial

4. The swab must be placed in the circular hole in the base. When carrying out this step it is **VERY IMPORTANT** not to exert any pressure on the swab as it is inserted. Simply drop it in the hole without pressing (Figure 8a).

5. Press down the plastic buffer vial to break the buffer vial seal and then squeeze laterally the vial twice to release completely the liquid. (Figure 8b)

6. Once these steps have been carried out the Oral Fluid DrugSIP® Cassette sample is ready to be analysed by the DrugSIP® Multi Analyzer device. **The QR code must be visible before placing the kit inside the device.** The detection kit should look as shown in Figure 9. If the QR code is not visible, the screen shows the message “Check card” and the device must be restarted.


Figure 9. Position of the card to be inserted in the DrugSIP® Multi Analyzer

5.2.2 Processing Oral Fluid DrugSIP® Cassette

To process and read the test using the DrugSIP® Multi Analyzer device provided by AlphaSIP, follow these steps:


1. Switch on the DrugSIP® Multi Analyzer. If the temperature is lower than the optimum temperature to ensure the correct processing of the test, the display will show: **“Heating Wait...”**. This process takes around 10-15 minutes.

WARNING: The device allows to perform the test although the temperature is not optimum (20°C), but the device should only be used when the screen shows the message: “DrugSIP Ready for test”. It is recommended to switch on the DrugSIP® Multi Analyzer before using it if the room temperature is lower than 20°C.

2. Once the device is switched on and the system shows the **“DrugSIP Ready for test”** screen, the user can process and read the test.

3. The user must be inserted the Oral Fluid DrugSIP® Cassette with the swab containing the saliva sample and the buffer vial once it has been pressing down and it has been squeezed into the DrugSIP® Multi Analyzer device. Perform the following steps:

Slide the tab in the left-hand side of the device and insert the Oral Fluid DrugSIP® Cassette following the marks located on the internal platform. Once the Oral Fluid DrugSIP® Cassette is positioned correctly you must released the tab, allowing the kit to be positioned correctly. Figure 10 shows the sequence to properly position the Oral Fluid DrugSIP® Cassette. Then close the door of the device.



Figure 10. Correct disposition of the Oral Fluid DrugSIP® Cassette in the DrugSIP® Multi Analyzer

WARNING: Do not run the test with the door open to avoid errors.

4. Then press the **BLUE BUTTON** to start the test. The screen will show “Identifying test” and after a few seconds “Oral Drug Test identified” will appear.

Later, the message “Simple-> blue / Multi-> yellow” will show in the screen.

In this case explained a simple test so that the **BLUE BUTTON** have to be pushed.

5. Then press the blue button so that the processing of the test can start. The message “**STEP 1 Loading ...**” will appear on the screen.

The first step takes around 3 minutes. When this step finishes, the following message will appear: “**STEP 1 Finished**” and then the device processes the saliva sample.


6. Automatically the second processing will start. Do not press any button. A several minutes countdown will begin. At the end of this countdown the processing of results will start. During this process the screen will show “**Analyzing...**”.

WARNING: If the device does not provide the results and the message “Analyzing...” does not disappear from the screen the equipment must be switched off, switched on again and the yellow button must be pressed.

If necessary, the device will carry out a second processing. The screen will show the message “**Double test**” and later “**Verifying...**”.


If after the verification “**Invalid test**” appears on the screen again, the device will print a read-out if no action is carried out. This will show that definitive results of the test have not been obtained.

WARNING: In the event that “Invalid test” appears, put the test back in the DrugSIP® Multi Analyzer and carry out the reading of processed tests (Point 9) to verify that the results are not valid.

7. When the analysis process has finished the message “**Extract test**” will appear on the screen. The user can now open the door and remove the test.

8. The results obtained will be displayed. The results are displayed on the screen alternately and repeatedly by means of the initials of the analysed substances, accompanied by the symbols “+” or “-” according to whether the results indicate the presence or absence of drugs, respectively.

For further information consult section 5.4 Coding of results.

The results will be printed out once. While the device is printing the screen shows the message “**Printing...**”

9. The test number will appear on the results print-out. This number must be written on the Oral Fluid DrugSIP® Cassette underneath the name of the person who has taken the test. The printed ticket will show the following information:

- Device information: device model and serial number.

- Test information: Date, time, test number, place to write subject’s identification, test interpretation indicating the initials of analysed drugs, the result (positive or negative) and finally another area for the subject to sign.



Figure 11. Oral Fluid DrugSIP® Cassette ticket

10. When the printing of results is complete the device will go back to default mode. The message “**DrugSIP Ready for test**” will appear on the screen and the test can be considered complete.

The results of the tests are stored in the internal memory of the DrugSIP® Multi Analyzer and they can be consulted using a PC and the data processing software included in the carrying case.

WARNING: If the device does not complete the whole reading cycle, test its battery level.

WARNING: The device is set up to read only the official Oral Fluid DrugSIP® Cassette cards provided by AlphaSIP.

5.3 Performing a simple urine test

5.3.1 Collecting an urine sample

To collect the urine sample, use the DrugSIP® urine collection kit (Figure 12b), which consists of: sterile container, thermometric strip and protection sticker. Once the sample has been collected, use the DrugSIP® urine detection kit (Figure 12a), which consists of: diagnosis card, base and plastic pipette. To perform the urine test correctly follow these steps:

Carefully read the instructions enclosed in the Urine DrugSIP® Cassette before use. They may include improvements not described in this manual


Figure 12. Urine DrugSIP® Cassette (a) DrugSIP® urine detection kit (b) DrugSIP® urine collection kit

1. The sample collection container should be opened in the presence of the person carrying out the test in order to verify its integrity. Remove the thermometric strip from the paper and position it on the lower part of the sample container, ensuring that it is attached well, as shown in Figure 13.

Figure 13. Sample collector container with the thermometric strip

2. Give the sample container to the subject in order to collect the urine. While the subject collects the sample, his or her identification must be written on the front of the diagnosis card in the space allocated for this (Figure 14).

Figure 14. Allocated space


Figure 15. Changes in colour spheres

3. Verify that the urine sample has not been altered or manipulated after collection. Check that changes in colour take place (blue-green) in the spheres on the thermometric strip, as shown in Figure 15. Otherwise, avoid using that sample for any test. (The sample can be stored at 4-8 °C if the test does not be performed in that moment.)

4. Once the validity of the sample is verified, there are two different methods to perform the test:

A. Immerse the diagnosis card in the container with the sample. Insert the card up to the line (DIP LINE) for 3 to 5 seconds. In figure 16a, you can see the correct height at which the card should be inserted in the sample.

B. Using the plastic pipette included in the kit, add 3 drops of urine (120ul) to the sample well located at the end of the diagnosis card. Figure 16b shows the correct place to add the sample.


Figure 16. Dispensing the sample
(a) Immerse the diagnostic card (b) Use the plastic pipette

WARNING: If the card is inserted in the sample beyond the test DIP LINE, or the sample is pipetted directly over the strips, the result of the test will not be valid.

5. Put the protection sticker on the bottom of the card to prevent the sample dripping (Figure 17a). Insert the diagnosis card in the plastic base provided, holding it by the plastic edges to avoid touching the reactive strips, as shown in Figure 17b.


Figure 17. Inserting the card in the plastic base
(a) Protection sticker to proof the contact area (b) Plastic base for reading the test

6. Once these steps have been carried out, the Urine DrugSIP® Cassette sample is ready to be analyzed by the DrugSIP® Multi Analyzer device. **The QR code must be visible before placing the kit inside the device.** The detection kit should be as shown in Figure 18. If the QR code is not visible, the screen shows the message “Check card” and the device must be restarted.


Figure 18. Correctly placing the urine sample card in the DrugSIP® Multi Analyzer

5.3.2 Processing Urine DrugSIP® Cassette

To process and read the test using the DrugSIP® Multi Analyzer device provided by AlphaSIP, follow these steps:

1. Switch on the DrugSIP® Multi Analyzer. If the temperature is lower than the optimum temperature to ensure the correct processing of the test, the display will show: **“Heating Wait...”**. This process takes around 10-15 minutes.

WARNING: The device allows to perform the test although the temperature is not optimum (20°C), but the device should only be used when the screen shows the message: “DrugSIP Ready for test”. It is recommended to switch on the DrugSIP® Multi Analyzer before use it if the room temperature is lower than 20°C.

2. Once the device is switched on and the system shows the **“DrugSIP Ready for test”** screen, the user can process and read the test.

3. The user must be inserted the Urine DrugSIP® Cassette with the diagnosis card containing the urine sample and protected with the sticker. The card must be inserted in the base and then into the DrugSIP® Multi Analyzer device. Do the following steps:

Slide the tab in the left-hand side of the device and insert the Urine DrugSIP® Cassette following the marks located on the internal platform. Once the Urine DrugSIP® Cassette is positioned correctly you must released the tab, allowing the kit to be positioned correctly. Figure 19 shows the sequence to properly position the Urine DrugSIP® Cassette. Then close the door of the device.

WARNING: Do not run the test with the door open to avoid errors.



Figure 19. Correctly placing the Urine DrugSIP® Cassette in the DrugSIP® Multi Analyzer

WARNING: Do not run the test with the door open to avoid errors.

4. Then press the **BLUE BUTTON** to start the test. The screen will show **"Identifying test"** and after a few seconds **"Urine Drug Test identified"** will appear.

5. Once identified, the device will start a countdown to process the test.

6. At the end of this countdown the processing of results will start. During this process the screen will show **"Analyzing..."**.

WARNING: If the device does not provide the results and the message "Analyzing..." does not disappear from the screen the equipment must be switched off, switched on again and the yellow button pressed.

If necessary, the device will carry out a second processing. The screen will show the message **"Double test"** and later **"Verifying..."**.


If after the verification **"Invalid test"** appears on the screen again, the device will print a read-out if no action is carried out. This will show that definitive results of the test have not been obtained.

WARNING: In the event that "Invalid test" appears, put the test back in the DrugSIP® Multi Analyzer and carry out the reading of processed tests (Point 9) to verify that the results are not valid.

7. When the analysis process has finished the message **"Extract test"** will appear on the screen. The user can now open the door and remove the test.

8. The results obtained will be displayed. The results are displayed on the screen alternately and repeatedly by means of the initials of the analysed substances, accompanied by the symbols "+" or "-" according to whether the results indicate the presence or absence of drugs, respectively.

For further information consult section 5.4 Coding of results.

The results will be printed out once. While the device is printing the screen shows the message “**Printing...**”

9. The test number will appear on the results print-out. This number must be written on the Urine Fluid DrugSIP® Cassette underneath the name of the person who has taken the test (Figure 20). The printed ticket will show the following information:

- Device information: device model and serial number.

- Test information: Date, time, test number, place to write subject's identification, test interpretation indicating the initials of analysed drugs, the result (positive or negative) and finally another area for the subject to sign.


Figure 20. Urine DrugSIP® Cassette ticket

10. When the printing of results is complete the device will go back to default mode. The message “**DrugSIP Ready for test**” will appear on the screen and the test can be considered complete.

The results of the tests are stored in the internal memory of the DrugSIP® Multi Analyzer and they can be consulted using a PC and the data processing software included in the carrying case.

WARNING: If the device does not complete the whole reading cycle, test its battery level.
WARNING: The device is set up to read only the official Urine DrugSIP® Cassette cards provided by AlphaSIP.

5.4 Coding of results

Once the device has finished the reading process, the results will be displayed on the screen of the device as follows.

Code	Name
AMP	Amphetamine
BAR	Barbiturates
BZO	Benzodiazepines
COC	Cocaine
MTD	Methadone
MET	Metamphetamine
OPI	Opiates
PCP	Phencyclidine
PPX	Propoxyphene
THC	Cannabis/Marijuana

In the second line, just underneath the code for each drug, the test results will be displayed as they appear in the following table:

Message on screen	Result
-	Negative (absence of drug in a higher concentration than cut off)
+	Positive (presence of drug in a higher concentration than cut off)
Invalid test	No conclusive test (See troubleshooting)

6. Temperature control of the sample

A temperature control system for the sample is integrated into the DrugSIP® Multi Analyzer device. This way the device guarantees that the test reading is carried out at an optimum temperature of 20°C, ensuring correct processing of the test, regardless of the outside temperature.

The device heats the area to be used for the DrugSIP® Cassette depending on the room temperature. It also adapts in real time to the reference temperature while processing the test.

Whenever the device is at a room temperature of more than 20°C, the system will not start working.

WARNING: The device allows to perform the test although the temperature is not optimum (20°C), but the device should only be used when the screen shows the message: “DrugSIP Ready for test”. It is recommended to switch on the DrugSIP® Multi Analyzer before using it if the room temperature is lower than 20°C.

7. Software PC

7.1 Introduction

DrugSIP® Multi Analyzer allows you to download the stored tests on a PC, so that you can keep the results for future reference, consult the data of the realized tests and obtain statistics.

The device includes an internal memory capable of storing up to 10.000 tests results. In order to access to these files, the device must be connected to the computer using the supplied USB cable.

The data downloading works both if the device is switched on or not. Once connected to the PC, the files can be consulted as if it was a conventional USB memory.

To use the data downloading feature, you must install the software application included on the provided CD.


7.2 Software installation

In order to run the software properly, you must have installed the “.NET framework 4.0” package or higher on your PC. If you already have installed it, omit this step. If not, you can download it for free from the official webpage of Microsoft. Besides, this package is included on the provided CD (file “dotNetFx40_Full_x86_x64.exe” in the folder “compl”).

NOTE: Before the installation, it is highly recommended to deactivate the antivirus software temporarily. Otherwise, the installation process cannot be warranted to be completed successfully.

When inserting the CD in the optical device of the computer, if the autorun feature is enable in the operative system, the installation program will be launched automatically. If this feature is not enabled, it must be done manually, double clicking on the “Setup.exe” located on the root of the CD.


When the installation software is launched the next screen is shown:

First, you must select the desired language for the software. You can choose between English or Spanish. Then, you must select the installation directory for the software to be installed. This directory can be a default one, which is offered by the installation program (C:\DrugSIP_v2.3), or a directory chosen by the user. Once the desired language and directory options are selected, the program will be installed clicking on “Install”.

IMPORTANT NOTE: if you choose a directory which requires administrator rights to accede, once the program is installed you will have to execute it as administrator. In case of doubt, and whenever possible, it is recommended to use the default directory for the installation.

When the installation ends, a desktop shortcut is generated, which gives access to the data downloading software. In this moment you can turn on the antivirus software again. Depending on the operative system, the next screen could be shown. If this occurs, simply click on “This program has been correctly installed”.


7.3 Application use

The application is executed double clicking on the desktop icon, or double clicking the “DrugSIP.exe” file in the previously selected installation directory.

When the application starts, the next options are available:

- Start download from DrugSIP Analyzer
- Analyze previously downloaded data
- Exit

Clicking on the “Start download from DrugSIP Analyzer” button, the communication between the DrugSIP® Multi Analyzer and the PC starts, and the stored tests are copied from the device to the PC. In order to do that, you must follow this sequence:


- Connect the DrugSIP® Multi Analyzer device to the PC using the USB cable.
- Wait until the system recognizes the device. The device will be shown on the PC as if it was a conventional data storage unit.
- You can click now on “Start download from DrugSIP Analyzer”. When the data downloading process ends, the program will automatically accede to the data treatment menu.

Pressing on “Analyze previously downloaded data”, you accede directly to the data treatment menu, using previously downloaded tests, if there is any. This menu allows to analyze and generate statistics from the existing data.

Finally, clicking on “Exit” the application closes.

7.3.1 Data analysis


After clicking on “Start download from DrugSIP Analyzer” or “Analyze previously downloaded data”, the data treatment menu is shown. The tests are divided in two columns: saliva tests and urine tests. If there was no saliva test, that column would not be displayed, and the same would happen for the urine tests.

To show the information from a test, you can select it clicking on its name on the list or on its checkbox. Subsequently, you must click on “Analyze selected test”. If you have a very long list of tests and you want to limit it, you can do it with the “Filter available tests by date” option and clicking on “Show tests”.


The results of the test are opened in a new window, where the following information is shown:

- On left side of the window there are two graphics. This graphics correspond to the lateral flow results of the reactive strips. The one marked as “Lateral Flow 1” corresponds to the right strip, and the one marked as “Lateral Flow 2” corresponds to the left strip.
- The photograph of the reactive strips is located on the central part of the window
- On the right side there are the results of the test, as well as the date and time, the test number and the serial number of the DrugSIP® Multi Analyzer. If the test is valid, the result for each drug will be either “Positive” or “Negative”. If the test is invalid, the result will be “Null” for each drug.

Apart from the previous data, the window shows three buttons. The first one is “Generate report”, the second one is “Return” (which is used to return to the data treatment menu) and the third one is “Exit” (which closes the application).

Clicking on “Generate report”, a new window shows, where the user can fill in some data of the test, as shown below:


When clicking on “Continue”, the program elaborates a report in PDF format, including the information of the test (data, results, photography and graphics), in the desired directory of the PC.

7.3.2 Statistics

The statistics option allows the user to obtain a statistical summary of the tests. Clicking on “Saliva statistics”, the statistics corresponding to saliva tests will be displayed, and the same for “Urine statistics”.

When clicking on the corresponding button, a new window appears, where a data range can be established, as well as selecting all the available tests.

The statistics which are related to the selected tests are displayed pressing on “Show statistics”:

8. Troubleshooting

Numbering of errors	Potential Issues	Solutions
Error 01: Mechanism (hardware). Referring to mechanical pieces which push down the different components of the kit (buffer cartridge, swab, diagnostic card). In order to identify Error 01, you must open the device door during STEP 1.	Buffer cartridge actuator locking *Actuator: Mechanical part which pushes down some component of the kit towards the base	Please, shut down the device and turn it on again. When the actuators are at the top, remove the DrugSIP® Cassette. Put it again on the system platform ensuring that it is in the correct position. The base of the cassette must be perfectly aligned with the base guides. Start the test. If this solution does not work, please contact our technical service.
	Oral swab actuator locking	
	Diagnostic card actuator locking	
Check card (Restart device) The device has not been able to identify the type of test inserted	The QR code is not visible	Restart the device pressing the on/off button.
	The position of the QR code is incorrect	Ensure that the QR code is not damaged. Insert the DrugSIP® Cassette in the correct position, with the QR code visible. Insert the test in the device again and start the test.
	The QR code is damaged	If the QR code is damaged, substitute the test
Calibration done (Restart device)	The device has performed an internal auto-calibration procedure and it must be restarted manually.	Shut down the device and turn it on again to start the system initialization. If this solution does not work, please contact our technical service.
Non printing The Drugsip cannot print results on the ticket. (Frequently asked question num.5 - pg 43)	Problems with the paper of the ticket	If the printer's green led is blinking. Check the thermal paper roll or replace it in the event that it is finished.
	Hardware problems	If the printer's green led is turned off. Check the device's battery level and charge the device if it is necessary. With the red indicator led turned on, press DrugSIP's yellow button to analyze the test again. At the end of the analyzing process, the device must print the ticket. If this solution does not work, please contact our technical service.
Data downloading PC-DrugSIP connection	Problems downloading data	A) Verify the USB connection and try to connect with the PC again B) Reinstall Data Downloading program If these solutions do not work, please contact our technical service.
Analyzing (Block)	The device does not finish the analysis step	Restart the device pressing the on/off button. Press the yellow button for reading the test.
Invalid Test The result shown on the screen is NCL (The control lines of the strips are not detected or read)	Immune-Reactive strips do not develop: Problems with buffer solution.	Please, make sure that the Cassette is in the required position and the buffer vial has been checked and released previously. Then press the yellow button to analyze the test. If the “Invalid Test” error persists, you must repeat the test with a new cassette
	Immune-Reactive strips with slow developing process	
	Immune-Reactive developed	
	Moved or displaced Cassette	

9. Reading developed tests

The DrugSIP® Multi Analyzer device includes the option to read tests which have been carried out previously. To only read a test follow these steps:

WARNING: If more than 12 minutes have passed between the end of the incubation period and the start of the reading process (Oral Fluid DrugSIP® Cassette), the results obtained are not reliable.

1. Place the Oral Fluid DrugSIP® Cassette or Urine DrugSIP® Cassette you would like to read inside the device. The card should be entered as shown in Figure 21. **The QR code should be visible when placing the kit inside the device.**

Oral Fluid DrugSIP® Cassette

Urine DrugSIP® Cassette

Figure 21. Oral Fluid DrugSIP® Cassette and Urine DrugSIP® Cassette ready to be read.

2. Close the door and press the **YELLOW BUTTON**. The screen will show **"Identifying..."** and after a few seconds **"Oral Drug Test identified"** or **"Urine Drug Test identified"** will appear depending on the kind of test that the device is reading.

3. The device will read the test. The message **"Analyzing..."** will appear on the screen.

WARNING: If the device does not provide the results and the message "Analyzing..." does not disappear from the screen the equipment must be switched off, switched on again and the yellow button pressed.

If necessary, the device will carry out a second processing. The screen will show the message **"Double test"** and later **"Verifying..."**.

If after the verification **"Invalid test"** appears on the screen again, the device will print a read-out if no action is carried out. This will show that definitive results of the test have not been obtained.

WARNING: In the event that "Invalid test" appears, put the test back in the DrugSIP® Multi Analyzer and carry out the reading of processed tests (Point 9) to verify that the results are not valid.

4. When the analysis process has finished the message **"Extract test"** will appear on the screen. The user can now open the door and remove the test.

5. The results obtained will be displayed. The results are displayed on the screen alternately and repeatedly by means of the initials of the analysed substances, accompanied by the symbols "+" or "-" according to whether the results indicate the presence or absence of drugs, respectively.

The results will be printed out once. While the device is printing the screen shows the message **"Printing..."**

6. When the printing of results is complete the device will go back to default mode. The message **"DrugSIP Ready for test"** will appear on the screen and the test can be considered complete.

10. Safety

- The DrugSIP® Multi Analyzer drug detection device must be used by previously trained personnel who know the equipment and how to handle it according to the user manual.
- Establish all measures required by current regulations regarding safety and hygiene at work.
- When using the device the user must choose the personal protection equipment which corresponds to the hazard category presented by the fluid to be handled. If appropriate equipment is not used or is faulty the user will be exposed to liquids which could splash or spill.
- Do not place the DrugSIP® Multi Analyzer drug detection device in areas near to sources of heat (lighters, radiators, etc).
- Ensure that the door of the device is completely closed before performing the test so that no body parts (hair, clothes or jewellery) may get caught in moving parts.
- For any handling, cleaning, checking or substitution of any internal part of the device contact the manufacturer.
- Do not attempt to repair the device yourself; apart from invalidating the guarantee, this could also lead to damage in the general functioning of the equipment, as well as personal injuries (burns, wounds, damage, etc.) and damage to the electrical system.
- The device must never come in contact with water. If for any reason you suspect that water or any liquid has entered the device switch it off immediately and contact the manufacturer.
- Manufactured in accordance with European Directives for electrical safety, electromagnetic compatibility and machine safety.
- The device is recommended for use on stable surfaces to prevent falls which may damage device.

11. Maintenance and storage

In order for the detection kit to function properly it is necessary to follow certain recommendations.

NOTE: All the above mentioned rules for use will have no value if the device is not maintained regularly.

- Follow the instructions and warnings relating to this manual.
- Always keep this manual at hand so that anyone can consult it.
- Always use original components and replacements parts. Other devices may look similar but using them could damage the equipment.
- In the event of failure contact the manufacturer. Do not attempt to repair the device yourself, as apart from invalidating the guarantee, this could lead to damage in the general functioning of the equipment.
- For proper preservation of kit components it is recommended that they are kept at a temperature of between 10 and 35 °C, in a cool atmosphere and not exposed to sunlight. It is strongly recommended that all kit components are kept away from any heat source.

12. Cleaning

(Waste electrical and electronic equipment) (Applicable in the European Union and in European countries with systems of selective collection of waste.)

The presence of this symbol on the product, accessories or information material accompanying it indicates that at the end of its useful life, neither the product nor its electronic accessories (charger, cables, etc.) should be disposed of with other household waste. To avoid possible damage to the environment or human health represented by uncontrolled waste disposal, separate these products from other types of waste and recycle them correctly. This promotes sustainable reuse of material resources.

Individual users can contact the establishment where they bought the product or the appropriate local authorities to find out how and where to leave it so that it can be recycled safely and in an environmentally friendly way.

Commercial users can contact their supplier and consult the conditions of the purchase contract. This product and its electronic accessories should not be disposed of with other commercial waste.

Frequently asked questions (FAQ)

1- What should I do if the device does not turn on?

In the event that the device does not turn on after activating the ON/OFF switch, connect the battery charger to the device and wait for 30 minutes. After this time has elapsed, if the device still does not turn on, contact the technical service of AlphaSIP.

2- What should I do if the device does not complete the analysis or verification process?

If the device is blocked when the screen shows "Analyzing..." or "Verifying..." and the progress charge bar has not moved forward for longer than 2 minutes, shut down the device and check the battery level. If the red indicator light comes on, connect the device to a power source. In the event that this action does not solve the problem, contact the technical service of AlphaSIP.

3- What should I do if the device does not show the results during it is connected to the PC?

Shut down the device and disconnect from the PC. Turn on the device and press the appropriate button for reading the test. If the screen does not show the results is due to the device does not allow to perform a test while it is connected to the PC. In the event that this action does not solve the problem, contact the technical service of AlphaSIP.

4- When should I charge the battery of the device?

The device has a low battery indicator (red indicator light). If it comes on, connect the device to a power source.

In the case of an infrequent use of the device, it is recommended to charge it once a month to make sure a good battery maintenance.

5- Why does the device not print the final ticket?

First, check if during the printing process the green indicator light of the printer is blinking or it remains off. If the green indicator light remains off, check the battery level and connect the device to a power source if necessary. In the event that the green indicator light blinks, check if the device has printer paper roll. If this action does not solve the problem, contact the technical service of AlphaSIP.

6- What should I do if the printer paper becomes jammed?

In the event that the paper gets jammed in the printer, open the printer door by gently pulling the lever just below the led and remove the piece of paper that is jammed (Point 4.2.2). If the printer still does not work after removing the paper, contact the technical service of AlphaSIP.

7- How to change language, date and time?

Turn on the DrugSIP® Multi Analyzer while holding down the blue button for at least 5 seconds until the screen shows **“CHANGE LANGUAGE DATE AND TIME”** (Point 4.2.1).

8- What should I do if date and time do not refresh?

If the device's date and time do not refresh, contact the technical service of AlphaSIP.

9- The obtained results show NCL. ¿What should I do?

Verify that the control line has appeared at the top of the strips. If so, contact the technical service of AlphaSIP.

10- The screen shows the message “Run test configuration”. ¿What should I do?

If the screen shows the previous message or the device is in the configuration mode, please shut down the device and contact the technical service of AlphaSIP.

11- What should I do if data downloading software does not recognize the device?

The PC may need a few seconds to check up the new system. Please, wait for a moment and restart data downloading process. If the PC does not recognize the device, reconnect the USB cable and observe if some blue indicator lights appear near the cable connection. In the event that this action does not solve the problem, contact the technical service of AlphaSIP.

12- What should I do if the data downloading software shows an error?

If the PC shows a pop-up window with an error message during data downloading process, contact the technical service of AlphaSIP.

13- Why does the green indicator led not turn on when the device is connected to a power source?

The device may have suffered a battery over-discharge and therefore it needs a long charging period to achieve an adequate performance level. Charge the device for several hours and visualize if the green indicator light turns on. In the event that the green indicator light turns on, keep the device charging until it reaches the maximum level. On the contrary, if the green indicator light does not turn on, please contact the technical service of AlphaSIP.

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©2014 All rights reserved. AlphaSIP, DrugSIP and their logos are registered brands from Laboratorios Alpha San Ignacio Pharma. DSIP - D0315