

GENEDIA Quantum FR101

GFR-100



GENEDIA Quantum FR101 uses the GENEDIA Quantum test kit to measure results through in vitro diagnostic disposable tests and quantitative and qualitative fluorescence immunoassay.

HIGH ACCURACY



COMPACT SIZE



FAST RESULT



MULTI-USE



GENEDIA Quantum

FR101



GFR-100

- Provide objective results with an automated interpretation system
- Compact size optimized for on-site diagnosis
- Automatically recognize Operator ID, Patient ID, and Specimen ID using a built-in barcode scanner
- Report and store results in multiple ways to meet user's data management needs

SPECIFICATION

Product name	GENEDIA Quantum FR101	Power	"100-240V AC – 0.58A." for the AC/DC Adapter "5V DC – 4.0A" for the equipment
Measurement method	Fluorescence immunoassay	Display	4.3 inch LCD with touch screen
Running time	Test Mode 20 minutes, Scan Mode 20 seconds	Operating / Storage temperature	17 ~ 32°C / -10 ~ 60°C
Memory	1,000 Test	Operating / Storage humidity	15~75%
Weight	1,000g		
Dimensions	117(W) x 144(L) x 138(H) mm		



GREEN CROSS MEDICAL SCIENCE
 26, Mugeuk-ro 65beon-gil, Geumwang-eup, Eumseong-gun,
 Chungcheongbuk-do, 27632 Republic of Korea
gcms.overseas@gccorp.com
www.greencrossms.com

MADE IN KOREA